

Flexible Endoscopy Suction Pump

VACO^{SP120}



The Suction pump **VACO^{SP120}** is dedicated to endoscopy procedures. It is designed with a special system of noise and vibration dampening for an extremely quiet and smooth operation with an incredibly low noise level of 39,4 dB. The pump is a maintenance-free, reliable and vibration-free device.

The main case is made of very sturdy, unbreakable, fire retardant material and the pump body itself is made of aluminum alloy. The device has a cooling system with extra fan and the color coded gauge enables a safe vacuum management. Furthermore, the pump is equipped with overflow protection and safety vacuum regulator.

- Continuous run
- Super silent, only 39.4 dB
- Vibration free
- Over flow protection

VACO^{SP120} EndoMed Systems Suction Pump

Item No.	Description
643120	Suction Pump, VACO SP120

Flexible Endoscopy Suction Pump



VACO SP120

Technical Specifications

Device	VACO SP120
Dimensions (WxHxD)	440 x 300 x 310 mm
Weight	12 Kg
Power Supply	230V / 50-60 Hz
Safety Class	Protection class: I Application part: Type BF
CE Standard	Classification acc. 93/42/EEC: IIa EN ISO 13485 93/42/EEC IEC 601-1, EN 60601-1C
Suction Capacity	38 l/min ±5%
Vacuum	0-90 kPa
Drainage	0-40 cmH2O
Noise Level	39.4 dB
Delivery contains	1x 2-liter suction bottle with holder, 2x sterile filters, 1x silicone tube set, 1x Power cable, user manual

Environmental Conditions

Temperature (Operation)	+10°C to +40°C
Temperature (Transport and Storage)	-20° C to +45°C
Max. relative humidity	During operation: 70% During storage and transport: 70%

Optional

Item No.	Description
643221	Suction Bottle 2 Liter
643222	Suction Bottle 2 Liter, Ring Holder
643223	Suction Bottle 2 Liter, Clamp
643224	Silicone Tube, 1 Meter
643225	Silicone Tube, 3 Meter
643210	Trolley for VACO SP120
643250	Filter for VACO SP120 (Qty. 100)